


URIGHT

PLUS Mini

BLOOD GLUCOSE MONITORING SYSTEM



$\pm 15\text{mg/dL}$ & $\pm 15\%$



- Pocket size
- Before / After meal (AC/PC) presettable
- Mini USB port
- Fast result with small blood sample

PLUS
Series

Clinical Data for URiGHT PLUS Mini



Summary

- The clinical accuracy study followed the guideline of ISO 15197: 2013 and EN ISO 15197: 2015.
- Clinical data demonstrated PLUS Mini met the current accuracy standards of ISO 15197: 2013 and EN ISO 15197: 2015.
- Out of total 960 tests, average 98.6% (> 95%) fell in the range of $\pm 15\text{mg/dL}$ or $\pm 15\%$ ($\pm 0.83\text{ mmol/L}$ or $\pm 15\%$) from the reference value.
- The study showed the accuracy performance of PLUS Mini monitoring system could be used on healthcare professionals and self-testing.

Capillary & venous test results

- When the sample < 100mg/dL, 100% of the sample fell within $\pm 15\text{ mg/dL}$
- When the sample $\geq 100\text{mg/dL}$, 98% of the sample fell within $\pm 15\%$

Test strip lot	n	$\pm 15\text{ mg/dL}$ or $\pm 15\%$ ($\pm 0.83\text{mmol/L}$ or $\pm 15\%$)	n	Within $\pm 15\text{ mg/dL}$ (Within $\pm 0.83\text{mmol/L}$)	n	Within $\pm 15\%$
Capillary blood testing						
Total	960	947 (98.6%)	312	312 (100%)	648	635 (98.0%)
TD13B105-B0I	320	316 (98.8%)	104	104 (100%)	216	216 (100%)
TD13D215-B0E	320	316 (98.8%)	104	104 (100%)	216	212 (98.1%)
TD13E107-B0D	320	315 (98.4%)	104	104 (100%)	216	211 (97.7%)
Venous blood testing						
Total	960	958 (99.8%)	312	302 (100%)	648	646 (99.7%)
TD13B105-B0I	320	320 (100%)	104	104 (100%)	216	216 (100%)
TD13D215-B0E	320	319 (98.7%)	104	104 (100%)	216	215 (99.5%)
TD13E107-B0D	320	319 (99.7%)	104	104 (100%)	216	215 (99.5%)



Results of error grid analysis demonstrated that 100% of the results fall in zone A which indicates the meters are clinical accurate.